

INSTRUCTIONS FOR MOUTHGUARD

Field of application

Playsafe triple is suitable for all sports. Three-layered sports mouthguard, hard intermediate layer.

Instructions for cleaning and maintenance

Application:

The material is slightly flexible but there is a hard layer between the soft layers. It means that insertion of placing is important and must follow the instructions of dentist in this.

After use:

- Wash well with water.
- Best is to thoroughly clean the inner and exterior side of the mouthguard with a toothbrush and soap. **Do not use any strongly perfumed soaps. Not suited are: tooth-paste (contains abrasive particles), mouth-wash (can cause discolouring) and water that is hotter than 50°C (deformation). Solvent-based cleaning agents cause delamination of multi-layered splints.**
- Again, wash well with water.
- Shake off the water or dry with a towel.
- Never blow-dry - deformation!
- Very important: Allow the mouthguard to completely dry! Keep at a dry place, at best in a box like that has aeration holes.
- Again wash with water before using it.

Generation of bad smells

If after some time the splint has adopted a strong smell, additionally put the splint for one hour in a non-perfumed, concentrated soap solution, afterwards thoroughly wash with water. Such soap suds remove most of the smell generating bacteria.

Discolouring

Soft thermoforming materials have the tendency to discolour. This intake of colour pigments can be reduced or avoided by careful maintenance, but it cannot be reversed. Mouthwashes and amalgam fillings can also cause discolouring.

Biocompatibility

All materials are tested on their biocompatibility and are physiologically harmless in their intended use. They are CE marked based on the EU regulation 2017/745 concerning medical products.